



E8-1155 Appleby Line, Burlington, ON L7L 5H9  
 Tel: (905) 335-7029 Fax: (905) 335-4154  
 Email: info@premierorthoticslab.com  
 www.premierorthoticslab.com

**Rush Order**  
**4 to 5 Days (\$25)**

**Separate Invoice**

OFFICE USE ONLY

Date \_\_\_\_\_

EA # \_\_\_\_\_

**CLINIC INFORMATION**

Clinic Name: \_\_\_\_\_  
 Practitioner Name: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 Phone: \_\_\_\_\_

**PATIENT INFORMATION**

Patient Name: \_\_\_\_\_  
 Age: \_\_\_\_\_  
 Weight: \_\_\_\_\_  
 Shoe Size: \_\_\_\_\_  
 Gender:  M  F

**TOP COVER OPTIONS**

- Ultra Hyde
- Vinyl
- Cushioned Cotton (ETC)
- Microcell Puff  blk  pink swirl  blk swirl  purple swirl  blue swirl  yellow swirl
- Soft Leather  Diabetic Top  Spenco (Add \$10)
- Sulcus
- Full
- 3/4"
- 1/8"
- 1/16"

\*Default Sport Topcover is 1/8" Cushioned Cotton


\*Default Dress Topcover is 1/16" Ultrahyde

\*Default Standard and Polypro Topcover is 1/8" Ultrahyde

**ORTHOTIC MODIFICATIONS**

- Metatarsal Pad  L  R  Full  1/8"  1/16"
- Heel Spur Pad  L  R
- Heel Cushion  L  R  1 mm  2 mm  3 mm
- Heel Lift  L  R  Removable  Attached
- Morton's Ext  L  R  soft  firm  rigid
- Reverse Morton's Ext  L  R  2 mm  3 mm
- Arch Fill  L  R  soft  med  hard
- Kinetic Wedge  L  R  2 mm  3 mm
- Neutral Heel Post  L  R  Varus  Left \_\_\_°  Right \_\_\_°  Valgus  Left \_\_\_°  Right \_\_\_°  firm  soft  Dual density
- Fore Foot Post  Medial  Left \_\_\_mm  Right \_\_\_mm  Lateral  Left \_\_\_mm  Right \_\_\_mm
- 1st Ray Cutout  L  R
- 1st Met Cutout  L  R
- Heel Cup  mild  moderate  deep
- Heel Hole  L  R
- PMP Pad  L  R  2 mm  3 mm

**ADDITIONAL INSTRUCTIONS**



**ORTHOTIC TYPE**

- STANDARD**  2 mm  3 mm  
 1/8" Ultrahyde topcover; Standard flex shell with 100% memory, vinyl underlay.
- POLYPRO**  2 mm  3 mm  
 Polypropylene shell for extra control; vinyl underlay; ideal for plantar fasciitis cases
- SPORT**  1 mm  2 mm  3 mm  
 Sport flex shell with red proform arch fill for maximum shock absorption and vinyl underlay.
- DRESS**  1 mm  2 mm  3 mm  
 Dress flex shell with vinyl underlay; Mild to moderate control low profile orthotic

- ACCOMMODATIVE**   
 Dual laminated EVA shell with soft arch fill, full length diabetic topcover, vinyl underlay. Extra depth footwear recommended.
- PREMIER ELITE**  3 mm  
 3mm Polypropylene shell with med arch fill, deep heel cup, vinyl underlay. Microcell Puff topcover unless otherwise indicated. This semi rigid podiatric orthotic is an aggressive motion control device.

- FULL EVA**  EVA  
 2 layer dual laminated EVA shell with full length 1/8" ultrahyde topcover, vinyl underlay.
- SLIM SOFT**   
 2 mm RCH shell with 1/8" poron arch fill; low heel cup; full diabetic topcover

- SANDAL/CLOG ORTHOTIC**   
 2 mm RCH shell; one piece 35 dur EVA liner with shallow heel; low arch profile  
 black top  brown top  grey top

- UCBL**   
 3 mm coloured RCH shell; deep heel; med/lat flanges and coloured puff top.

- GAIT PLATE**  
 Med/Lat flanges & deep heel cup



- SLIM DRESS**  1 mm  2 mm  
 Slim narrow grind shell, top cover with 1mm nyplex underlay, full length.  
 For non-removable liner shoes.

\*If no selection is made, default is Standard with 1/8" Ultrahyde and 1mm underlay.  
**Total # of Orthotics Ordered: \_\_\_\_\_ Total # of Shoes Ordered: \_\_\_\_\_**

**FOOTWEAR ORDER**

<b>1</b>	Make: _____
	Model: _____
	Size/Width: _____ <input type="checkbox"/> With Orthotic
<b>2</b>	Make: _____
	Model: _____
	Size/Width: _____ <input type="checkbox"/> With Orthotic
<b>3</b>	Make: _____
	Model: _____
	Size/Width: _____ <input type="checkbox"/> With Orthotic