



E8-1155 Appleby Line, Burlington, ON L7L 5H9
 Tel: (905) 335-7029 Fax: (905) 335-4154
 Email: info@premierorthoticslab.com
 www.premierorthoticslab.com

CLINIC INFORMATION

Clinic Name:
 Practitioner Name:
 Address:
 Phone:

PATIENT INFORMATION

Patient Name:
 Age: Weight:
 Shoe Size:
 Gender: M F

TOP COVER OPTIONS

- Ultra Hyde
- Vinyl
- Cushioned Cotton (ETC)
- Microcell Puff blk pink swirl blk swirl purple swirl blue swirl yellow swirl
- Soft Leather Diabetic Top Spenco (Add \$10)
- Sulcus
- Full
- 3/4"
- 1/8"
- 1/16"

*Default Sport Topcover is 1/8" Cushioned Cotton
 *Default Dress Topcover is 1/16" Ultrahyde
 *Default Standard and Polypro Topcover is 1/8" Ultrahyde

ORTHOTIC MODIFICATIONS

- Metatarsal Pad L R Full 1/8" 1/16"
- Heel Spur Pad L R
- Heel Cushion L R 1 mm 2 mm 3 mm
- Heel Lift L R Removable Attached
- Morton's Ext L R soft firm rigid
- Reverse Morton's Ext L R 2 mm 3 mm
- Arch Fill L R soft med hard
- Kinetic Wedge L R 2 mm 3 mm
- Neutral Heel Post L R Varus Left ___° Right ___° Valgus Left ___° Right ___° firm soft Dual density
- Fore Foot Post Medial Left ___mm Right ___mm Lateral Left ___mm Right ___mm
- 1st Ray Cutout L R
- 1st Met Cutout L R
- Heel Cup mild moderate deep
- Heel Hole L R
- PMP Pad L R 2 mm 3 mm

ADDITIONAL INSTRUCTIONS



Rush Order
4 to 5 Days (\$25)

Separate Invoice

OFFICE USE ONLY

Date _____

EA # _____

ORTHOTIC TYPE

- STANDARD** 2 mm 3 mm
 1/8" Ultrahyde topcover; Standard flex shell with 100% memory, vinyl underlay.
- POLYPRO** 2 mm 3 mm
 Polypropylene shell for extra control; vinyl underlay; ideal for plantar fasciitis cases
- SPORT** 1 mm 2 mm 3 mm
 Sport flex shell with red proform arch fill for maximum shock absorption and vinyl underlay.
- DRESS** 1 mm 2 mm 3 mm
 Dress flex shell with vinyl underlay; Mild to moderate control low profile orthotic

ACCOMMODATIVE
 Dual laminated EVA shell with soft arch fill, full length diabetic topcover, vinyl underlay. Extra depth footwear recommended.

PREMIER ELITE 3 mm
 3mm Polypropylene shell with med arch fill, deep heel cup, vinyl underlay. Microcell Puff topcover unless otherwise indicated. This semi rigid podiatric orthotic is an aggressive motion control device.

FULL EVA EVA
 2 layer dual laminated EVA shell with full length 1/8" ultrahyde topcover, vinyl underlay.

SLIM SOFT
 2 mm RCH shell with 1/8" poron arch fill; low heel cup; full diabetic topcover

SANDAL/CLOG ORTHOTIC
 2 mm RCH shell; one piece 35 dur EVA liner with shallow heel; low arch profile
 black top brown top grey top

UCBL
 3 mm coloured RCH shell; deep heel; med/lat flanges and coloured puff top.

GAIT PLATE
 Med/Lat flanges & deep heel cup



SLIM DRESS 1 mm 2 mm
 Slim narrow grind shell, top cover with 1mm nyplex underlay, full length.
 For non-removable liner shoes.

*If no selection is made, default is Standard with 1/8" Ultrahyde and 1mm underlay.

Total # of Orthotics Ordered: _____ Total # of Shoes Ordered: _____

FOOTWEAR ORDER

1 Make: _____
 Model: _____
 Size/Width: _____ With Orthotic

2 Make: _____
 Model: _____
 Size/Width: _____ With Orthotic

3 Make: _____
 Model: _____
 Size/Width: _____ With Orthotic