

**CLINIC INFORMATION**

Clinic Name: \_\_\_\_\_  
 Practitioner Name: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 Phone: \_\_\_\_\_

**PATIENT INFORMATION**

Patient Name: \_\_\_\_\_  
 Age: \_\_\_\_\_ Weight: \_\_\_\_\_ Shoe Size: \_\_\_\_\_  
 Gender:  M  F

**TOP COVER OPTIONS**

Ultra Hyde Only  1/8"  1/16" Profile:  Sulcus  
 Vinyl Only  1/8"  1/16"  Full  
 Cushioned Cotton (ETC)  3/4  
 Microcell Puff  blk  pink swirl  blk swirl  
 purple swirl  blue swirl  yellow swirl  1/8"  
 Soft Leather  Diabetic Top (Add \$10)  Spenco  1/16"

**ORTHOTIC MODIFICATIONS**

Metatarsal Pad  L  R  
 5mm  3mm  1.5mm

Heel Spur Pad  L  R

Heel Cushion  L  R  
 1 mm  2 mm  3 mm

Heel Lift  L  R  Removable  
 Attached  
 \_\_\_mm \_\_\_mm

Morton's Ext  L  R  
 To IPJ  End of Toes  soft  firm  rigid

Reverse Morton's Ext  L  R  
 2 mm  3 mm

Arch Fill  L  R  
 soft  med  hard

Arch Strip  L  R  
 Medial  Lateral  Hard  Soft

Kinetic Wedge  L  R  
 2 mm  3 mm

Neutral Heel Post  L  R  
 Extrinsic  Varus  Left \_\_\_°  Right \_\_\_°  
 Intrinsic  Valgus  Left \_\_\_°  Right \_\_\_°  
 firm  soft  Dual density

Fore Foot Post  Varus  Left \_\_\_°  Right \_\_\_°  
 Extrinsic  Valgus  Left \_\_\_°  Right \_\_\_°  
 Intrinsic

1st Ray Cutout  L  R

1st Met Cutout  L  R

Heel Cup  Flat  mild  moderate  deep

Heel Hole  L  R  
 1"  1.5"  2"

PMP Pad  L  R  Sulcus  
 2 mm  3 mm  End of Toes

Toe Crest  L  R

Flange/Flare  L  R  
 Medial  Lateral  Hard  Soft

Heel Skive  L  R  
 Medial  Lateral  2mm  4mm  6mm


Neuroma Pad  L  R  
 2-3  3-4  4-5

**GAIT PLATE**

Induce in toeing

Induce out toeing

**ADDITIONAL INSTRUCTIONS**

 \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_

**ADDITIONAL INSTRUCTIONS**

 \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_

**R L** \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_

**Separate Invoice** OFFICE USE ONLY  
 Date \_\_\_\_\_  
 **Rush Order** Foam  Plaster  Scan   
**4 to 5 Days (\$25)**

**ORTHOTIC TYPE**

**STANDARD**  1 mm  2 mm  3 mm  
 1/8" Ultrahyde topcover; Standard flex shell with 100% memory, vinyl underlay.

**POLYPRO**  2 mm  3 mm  
 Poly shell 1/8 ultra hyde top: vinyl underlay; ideal for plantar fasciitis cases

**SPORT**  1 mm  2 mm  3 mm  
 Sport flex shell with red proform arch fill 1/8 ETC top

**TURBO SPORTS**

Turbo Court  Turbo Soccer  Turbo Skate

Turbo Golf  Turbo Cycling  Turbo Run

**DRESS**  1 mm  2 mm  3 mm  
 Dress flex shell with vinyl underlay; 1/16 ultra hyde top

**ACCOMMODATIVE**   
 6mm EVA shell with soft arch fill, full length diabetic topcover, vinyl underlay. Extra depth footwear recommended.

**PREMIER ELITE**  3 mm  
 3mm Poly shell with med arch fill, deep heel cup, vinyl underlay. Microcell Puff topcover

**EVA**  3/4 LENGTH  FULL LENGTH EVA SHELL  
 9 mm EVA shell with full length 1/8" ultrahyde topcover, vinyl underlay.

**SLIM SOFT**   
 1 mm Poly shell with 1/8" poron arch fill; low heel cup; full diabetic topcover

**SANDAL/CLOG ORTHOTIC**   
 1 mm Poly shell; one piece 35 dur EVA liner with shallow heel; low arch profile, Full removable Liner  
 black top  brown top

**BIRKENSTOCK ORTHOTIC** (Non-exchangeable / returnable)  
 Built in Orthotic  Orthotic Separate (please specify)  
 Built in orthotic with EVA densities, No plastic shell

**UCBL**  1 mm  2 mm  3 mm  
 Deep heel; med/lat flanges and coloured puff top.

**SLIM DRESS**  1 mm  2 mm  3 mm  
 Slim narrow grind shell, suede top with 1mm nyplex underlay, full length.  
 For non-removable liner shoes.  
 \*If no flexion is made, default is Standard with 1/8" Ultrahyde and 1mm underlay.

**Total # of Orthotics Ordered: \_\_\_\_\_ Total # of Shoes Ordered: \_\_\_\_\_**

**FOOTWEAR ORDER**

**1** Make: \_\_\_\_\_  
 Model: \_\_\_\_\_  
 Size/Width: \_\_\_\_\_  With Orthotic

**2** Make: \_\_\_\_\_  
 Model: \_\_\_\_\_  
 Size/Width: \_\_\_\_\_  With Orthotic

**3** Make: \_\_\_\_\_  
 Model: \_\_\_\_\_  
 Size/Width: \_\_\_\_\_  With Orthotic

**COMPRESSION SOCKS ORDER**

Make	Model	Compression	Description	Sizes: S, M, L, XL	Qty:
Jiani	MK0058	20-30 mmHg	Black		
Jiani	MK0058	20-30 mmHg	Black / Ora		
Jiani	MK0058	20-30 mmHg	Navy / Blu		
Jiani	MK0058	20-30 mmHg	Blu / Grey		
Jiani	MK0058	20-30 mmHg	Blu / White / Yel		
FootNurse	JG - 360	20-30 mmHg	White		